





MAKAIRA®

INTERNAL
SPINAL
FIXATION
SYSTEM





MAKAIRA® FEATURES







The MAKAIRA® Toploading Pedicle Screw System is a low-profile and easy-to-use system that impresses the user with its simple design and good handling.

It is indicated for use in the thoracic and lumbar regions of the spine. The design of the screw with the large thread flanks and conical tip (a self-tapping thread) ensures that the screw can be placed very precisely and easily.

The range of screws and rods covers all common sizes for use in various indications. MAKAIRA® is a pedicle screw system, preferably used for open and small accesses.

Implants and instruments are supplied non-sterile packed and clearly placed on the set.


FEATURES:

-  4 Screw diameter
-  12 Screw lengths
-  Monoaxial and polyaxial screws
-  Polyaxial spondylosisthesis screws
-  Pre-bended rods in 11 lengths
-  Cross connectors



MAKAIRA® IMPLANTS

ARTICLE NUMBER STANDARD	DESCRIPTION	Ø	LENGHT	ARTICLE NUMBER SPONDY
700.303.845.0XX	Polyaxial Screw	4,2	25-80	
700.303.850.0XX	Polyaxial Screw	5,5	25-80	700.323.650.0XX
700.303.860.0XX	Polyaxial Screw	6,5	25-80	700.323.660.0XX
700.303.870.0XX	Polyaxial Screw	7,5	25-80	700.323.670.0XX
700.205.645.0XX	Monoaxial Screw	4,2	25-80	
700.205.650.0XX	Monoaxial Screw	5,5	25-80	
700.205.660.0XX	Monoaxial Screw	6,5	25-80	
700.206.670.0XX	Monoaxial Screw	7,5	25-80	



ARTICLE NUMBER	DESCRIPTION	Ø	Ø	LENGHT
700.101.550.XXX	Rod, Curved	5,5	-	25-125
700.101.550.XXX	Rod, Straight	5,5	-	130-400
700.403.700.106	Locking Screw			
700.403.700.102	Connecting Rod small			
700.403.700.104	Connecting Rod large			





MAKAIRA®

MORE INFORMATION

COMPATIBLE SYSTEMS

POSTERIOR OR TRANSFORAMINAL LUMBAR INTERBODY FUSION (PLIF CAGES) (TLIF CAGES)


- BASIC CAGE® PLIF
- PINA® Trochilus TLIF
- PINA® Expandus PLIF
- PINA® Plintus PLIF

Further Information at
info@bricon.com
and
+49 (0) 7461 9336-0

SIMILAR SYSTEMS

OCTOPODA®

- Self-breaking set-screw
- Ultra-Low profile of tulip
- More screw options
- Improved screw thread
- Increased screw angulation



BRICON

Bricon GmbH
Eisenbahnstrasse 100
78573 Wurmlingen
Germany

+49 7461 9336 - 0
info@bricon.com
bricon.com



[https://bricon.com/
produkte/wirbelsaeule/](https://bricon.com/produkte/wirbelsaeule/)

CE 0482