

Art.-Nr.	Width	Length	Height
CKK 12.14.04	12 mm	14 mm	4 mm
CKK 12.14.05	12 mm	14 mm	5 mm
CKK 12.14.06	12 mm	14 mm	6 mm
CKK 12.14.07	12 mm	14 mm	7 mm
CKK 12.14.08	12 mm	14 mm	8 mm
CKK 12.14.09	12 mm	14 mm	9 mm

Art.-Nr.	Width	Length	Height
CKK 14.16.04	14 mm	16 mm	4 mm
CKK 14.16.05	14 mm	16 mm	5 mm
CKK 14.16.06	14 mm	16 mm	6 mm
CKK 14.16.07	14 mm	16 mm	7 mm
CKK 14.16.08	14 mm	16 mm	8 mm
CKK 14.16.09	14 mm	16 mm	9 mm

Art.-Nr.	Width	Length	Height
CKK 12.12.05	12 mm	12 mm	5 mm
CKK 12.12.06	12 mm	12 mm	6 mm
CKK 12.12.07	12 mm	12 mm	7 mm

**sterile**

Art.-Nr.	Width	Length	Height
CKPOS 12.14.04	12 mm	14 mm	4 mm
CKPOS 12.14.05	12 mm	14 mm	5 mm
CKPOS 12.14.06	12 mm	14 mm	6 mm
CKPOS 12.14.07	12 mm	14 mm	7 mm
CKPOS 12.14.08	12 mm	14 mm	8 mm
CKPOS 12.14.09	12 mm	14 mm	9 mm

Art.-Nr.	Width	Length	Height
CKPOS 14.16.04	14 mm	16 mm	4 mm
CKPOS 14.16.05	14 mm	16 mm	5 mm
CKPOS 14.16.06	14 mm	16 mm	6 mm
CKPOS 14.16.07	14 mm	16 mm	7 mm
CKPOS 14.16.08	14 mm	16 mm	8 mm
CKPOS 14.16.09	14 mm	16 mm	9 mm

Art.-Nr.	Width	Length	Height
CKPOS 12.12.05	12 mm	12 mm	5 mm
CKPOS 12.12.06	12 mm	12 mm	6 mm
CKPOS 12.12.07	12 mm	12 mm	7 mm

# PINA®

*Your Swiss Option*

## CERVICUS *Invadur™*



### Description

The intervertebral cervical cages CERVICUS *Invadur*™ from PINA® Medizintechnik Vertriebs AG are implants intended for implantation between the end plates of the individual vertebrae. Their shape is adapted to the morphology of the intervertebral spaces being treated and the prescribed surgical techniques. The cages are available in various designs and sizes. The PINA® CERVICUS *Invadur*™ cages are intended for use in segments C2 to C7.



### Product characteristic

CERVICUS *Invadur*™ cages are made from \**Invadur*™. Due to the x-ray transparency there are no artefacts in MRI or CT. The x-ray markers made from Titanium allow the intraoperative and the postoperative control of the implant position. A large contact area allows a maximal bone contact and reduces the risk of subsidence. Through the existing openings, a filling with autologous bone or bone substitute material is possible. The CERVICUS *Invadur*™ is available sterile as well as non-sterile.



### Indications

- Degenerative disc diseases and instabilities
- Degenerative Spondylolisthesis
- Revision surgery

Changes reserved

\**Invadur*™ is the PINA-Brand for polymers out of the PEAK-group, such as PEEK and PEKK

4 degrees	Width	Length	Height
	12 mm	14 mm	4 mm
	12 mm	14 mm	5 mm
	12 mm	14 mm	6 mm
	12 mm	14 mm	7 mm
	12 mm	14 mm	8 mm
	12 mm	14 mm	9 mm

4 degrees	Width	Length	Height
	14 mm	16 mm	4 mm
	14 mm	16 mm	5 mm
	14 mm	16 mm	6 mm
	14 mm	16 mm	7 mm
	14 mm	16 mm	8 mm
	14 mm	16 mm	9 mm

4 degrees	Width	Length	Height
	12 mm	12 mm	5 mm
	12 mm	12 mm	6 mm
	12 mm	12 mm	7 mm

### Instruments CERVICUS *Invadur*™

