

PINA®

Your Swiss Option

CERVICUS


INVADUR®

ANTERIOR

CERVICAL

INTERBODY

FUSION






CERVICUS INVADUR

FEATURES

The intervertebral cervical CERVICUS Invadur cages from PINA® Medizintechnik Vertriebs AG are implants designed for implantation between the end plates of the vertebral bodies. These cages are intended for applications in segments C2 to C7. The implants can be inserted via a lateral anterior approach to the cervical spine. Their shape is adapted to the morphology of the intervertebral body spaces to be treated and to the prescribed surgical techniques. The cages are available in various designs and sizes. CERVICUS Invadur® cages are made of polymers from the PAEK group. This specially made Invadur® material is X-ray transparent, bio-, MRI-, and CT-compatible. CERVICUS Invadur cages are available both sterile and non-sterile.

FEATURES:

- ☐ Available in Sterile and non-sterile
- ☐ Lordotic angle of 4°
- ☐ Five different heights (4 mm - 9 mm) with an increment of 1 mm
- ☐ Easy-of-use instrumentation
- ☐ Five different Footprints
- ☐ Three-dot x-ray marker for visualization during imaging



CERVICUS INVADUR


IMPLANTS

NON-STERILE

ARTICEL NUMBER	ARTICLE DESCRIPTION	SIZE (WIDTH X LENGTH)	HEIGHT X	ANGLE
CKK 12.12.0X	CERVICUS Invadur	12 x 12	5 - 7	4°
CKK 12.14.0X	CERVICUS Invadur	12 x 14	4 - 9	4°
CKK 14.16.0X	CERVICUS Invadur	14 x 16	4 - 9	4°
CKK 17.14.0X	CERVICUS Invadur	17 x 14	5 - 7	4°
CKK 17.16.0X	CERVICUS Invadur	17 x 16	5 - 7	4°

STERILE

ARTICEL NUMBER	ARTICLE DESCRIPTION	SIZE (WIDTH X LENGTH)	HEIGHT X	ANGLE
CKPOS 12.12.0X	CERVICUS Invadur steril	12 x 12	5 - 7	4°
CKPOS 12.14.0X	CERVICUS Invadur steril	12 x 14	4 - 9	4°
CKPOS 14.16.0X	CERVICUS Invadur steril	14 x 16	4 - 9	4°
CKPOS 17.14.0X	CERVICUS Invadur steril	17 x 14	5 - 7	4°
CKPOS 17.16.0X	CERVICUS Invadur steril	17 x 16	5 - 7	4°





CERVICUS INVADUR

MORE INFORMATION

ALTERNATIVE SYSTEMS

→ Pina Echinus Invadur, Ti-cated




Further information at
info@bricon.com
and
+49 (0) 7461 9336-0

COMPATIBLE SYSTEMS

ANTERIOR CERVICAL PLATES

→ SHARK CP PRO® ACP

→ INION S1™ ACP



PINA®

Your Swiss Option

PINA Medizintechnik Vertriebs AG
Neuwiesenstrasse 15
8400 Winterthur / Switzerland

Tel.: +41 (44) 586 52 - 62

info@pina-med.ch
www.pina-med.ch

BRICON

Distributed by:

Bricon GmbH
Eisenbahnstrasse 100
78573 Wurmlingen / Germany

Tel.: +49 (0) 7461 9336 - 0

info@bricon.com
www.bricon.com

