

PINA[®]

Your Swiss Option

TROCHILUS


TITAN

TRANSFORAMINAL

LUMBAR

INTERBODY

FUSION











TROCHILUS TITAN

FEATURES

The PINA® TROCHILUS Titan TLIF-Cage System was developed for the unilateral supply of lumbar, degenerative changes after decompression and discectomy to secure the vertebral column, to restore intervertebral height and to secure lumbar lordosis. The broad, toothed surface profile provides good contact with the vertebral endplates whose integrity is to be maintained. For individual securing of the lumbar lordosis, different inclination angles were designed, plane-parallel and 8 degrees.

FEATURES:

-  Interlocking surface profile
-  Three different lengths (27 mm, 32 mm, 36 mm) and heights (8 mm, 10 mm, 12 mm)
-  Two different lordotic angles (0° und 8°)
-  Easy handling of the instrumentation
-  Wedge-shaped tip enables gentle implantation
-  Large filling openings for autogenous bonegraft



To fix the TROCHILUS cage firmly on the inserter for implantation, the nut is turned clockwise and tightened opposite the handle. Otherwise the cage is loosely, pivotally held on the inserter.



TROCHILUS TITAN

IMPLANTS

ARTICLE NUMBER	ARTICLE DESCRIPTION	IMPLANT LENGTH	IMPLANT HEIGHT	MAX. LORDOTIC ANGLE
TT 27.08.0.SL	TROCHILUS Titan	27	8	0°
TT 27.10.0.SL	TROCHILUS Titan	27	10	0°
TT 27.12.0.SL	TROCHILUS Titan	27	12	0°
TT 32.08.0.SL	TROCHILUS Titan	32	8	0°
TT 32.10.0.SL	TROCHILUS Titan	32	10	0°
TT 32.12.0.SL	TROCHILUS Titan	32	12	0°
TT 36.08.0.SL	TROCHILUS Titan	36	8	0°
TT 36.10.0.SL	TROCHILUS Titan	36	10	0°
TT 36.12.0.SL	TROCHILUS Titan	36	12	0°
TT 27.08.8.SL	TROCHILUS Titan	27	8	8°
TT 27.10.8.SL	TROCHILUS Titan	27	10	8°
TT 27.12.8.SL	TROCHILUS Titan	27	12	8°
TT 32.08.8.SL	TROCHILUS Titan	32	8	8°
TT 32.10.8.SL	TROCHILUS Titan	32	10	8°
TT 32.12.8.SL	TROCHILUS Titan	32	12	8°
TT 36.08.8.SL	TROCHILUS Titan	36	8	8°
TT 36.10.8.SL	TROCHILUS Titan	36	10	8°
TT 36.12.8.SL	TROCHILUS Titan	36	12	8°






TROCHILUS TITAN

MORE INFORMATION

COMPATIBLE SYSTEMS

BRICON OCTOPODA®

- Self-breaking set-screw
- Ultra-Low profile of tulip
- Multiple screw options
- Improved screw thread
- Increased screw angulation



Further information at

info@bricon.com

and

+49 (0) 7461 9336-0

PINA®

Your Swiss Option

PINA Medizintechnik Vertriebs AG

Neuwiesenstrasse 15

8400 Winterthur / Switzerland

Tel.: +41 (44) 586 52 - 62

info@pina-med.ch

www.pina-med.ch

CE 0297

BRICON

Distributed by:

Bricon GmbH

Eisenbahnstrasse 100

78573 Wurmlingen / Germany

Tel.: +49 (0) 7461 9336 - 0

info@bricon.com

www.bricon.com