




BRICON

SHARK CAGE®

A NTERIOR
C ERVICAL
I NTERBODY
F USION










SHARK CAGE®


FEATURES

SHARK CAGE®, the essential part of bricon cervical spine portfolio, used for anterior cervical interbody fusion. It covers extensive range of specification, including two shapes of footprints, four different combinations of width and depth, and eight heights, which could be powerful tool for surgeons to address the challenges posed by everyday surgery.

SHARK CAGE® is made of biocompatible PEEK material, can be used as stand-alone cage, or to be physically combined with the shark plate to form much stronger construct, or could be used with SHARK CP® or SHARK CP PRO®, two generations of ACP plate of Bricon.

FEATURES


-  Lordotic and Anatomic shape
-  Eight different heights
-  Teeth on the implant surface provide initial stability
-  Easy-of-use instrumentation
-  Three-point Radiopaque markers for visualisation during imaging




SHARK CAGE®

IMPLANTS

ARTICLE NUMBER	DESCRIPTION	SIZE(L x B)	HEIGHT (XX)
700.515.214.1XX	SHARK CAGE lordotic shape	13 x 14	5-12
700.515.217.1XX	SHARK CAGE lordotic shape	13 x 17	5-12
700.517.217.1XX	SHARK CAGE lordotic shape	16 x 17	5-12
700.517.219.1XX	SHARK CAGE lordotic shape	16 x 19	5-12



ARTICLE NUMBER	DESCRIPTION	SIZE (L x B)	HEIGHT (XX)
700.522.214.1XX	SHARK CAGE anatomic shape	13 x 14	5-12
700.522.217.1XX	SHARK CAGE anatomic shape	13 x 17	5-12
700.523.217.1XX	SHARK CAGE anatomic shape	16 x 17	5-12
700.523.219.1XX	SHARK CAGE anatomic shape	16 x 19	5-12



ARTICLE NUMBER	DESCRIPTION	Ø	LENGTH (XX)	HEIGHT (XX)
700.515.114.114	SHARK CAGE Plate			5-12
700.515.111.412	SHARK CAGE Screw	4	12-18	





SHARK CAGE®

MORE INFORMATION

INSTRUMENTS

- Simple insertion instrument
- Easy to use
- Self explaining (intuitiv)

Further information at
info@bricon.com
and
+49 (0) 7461 9336-0

COMPATIBLE SYSTEMS

ANTERIOR CERVICAL PLATES

- SHARK CP PRO® ACP
- SHARK CP® ACP
- INION S1™ ACP



BRICON

Bricon GmbH
Eisenbahnstrasse 100
78573 Wurmlingen
Germany

Phone +49 (0) 7461 9336-0

info@bricon.com
www.bricon.com

 0482

Distributed by: