

BRICON

SHARK CP PRO®

ANTERIOR
CERVICAL
PLATE














SHARK CP PRO[®]

FEATURES

A great leap on basis of the proven SHARK CP[®] system. SHARK CP PRO[®] has a completely new developed locking mechanism ensuring the optimum combination of safety and enhanced usability. The fully optimized instrumentation complement the implant system in setting a new clinical standard for anterior cervical plate.

FEATURES

-  Innovative locking mechanism
-  Slim waist of 14 mm, maximum width of 16 mm
-  41 Plate sizes, up to 5 Level, length ranging from 12 mm-105 mm
-  Ultra-low profile with 1,95 mm thickness
-  Screw angulation $\pm 10^\circ$, cranial/caudal up to 25°
-  Large windows for bonegraft application
-  Variable and fixed screws








SHARK CP PRO[®]

IMPLANTS

DESCRIPTION	ARTICLE NUMBER	Ø	LENGTH
Screw, variable, self-drilling	700.515.120.11X	Ø 4,0	10-18
Screw, variable, self-tapping	700.515.104.11X	Ø 4,0	10-18
Screw, variable, self-tapping	bi-cortical 700.515.104.12X	Ø 4,0	20-28
Rescue screw, var., self-tapping	700.515.116.41X	Ø 4,35	11-19
Rescue screw, var., self-tapping	bi-cortical 700.515.116.42X	Ø 4,35	21-29
Screwcaddy, variable screws	520.200.067	-	-
Lid for screwcaddy	520.200.069	-	-

DESCRIPTION	ARTICLE NUMBER	Ø	LENGTH
Screw, fix, self-drilling	700.515.115.41X	Ø 4,0	10-18
Screw, fix, self-tapping	700.515.126.11X	Ø 4,0	10-18
Screw, fix, self-tapping	bi-cortical 700.515.126.12X	Ø 4,0	20-28
Rescue screw, fix, self-tapping	700.515.113.41X	Ø 4,35	11-19
Rescue screw, fix, self-tapping	bi-cortical 700.515.113.42X	Ø 4,35	21-29
Screwcaddy, fixed screws	520.200.068	-	-
Lid for screwcaddy	520.200.069	-	-

DESCRIPTION	ARTICLE NUMBER	LENGTH	LEVEL
Anterior cervical plate 4	700.515.117.1XX	12-30mm	1 Level
Anterior cervical plate 6	700.515.117.1XX	30-44 mm	2 Level
Anterior cervical plate 8	700.515.117.1XX	42-66 mm	3 Level
Anterior cervical plate 10	700.515.117.1XX	64-88 mm	4 Level
Anterior cervical plate 12	700.515.117.151	75-105 mm	5 Level
Platecaddy	520.200.070	-	-
Lid for platecaddy	520.200.066	-	-






SHARK CP PRO[®]

MORE INFORMATION

INSTRUMENTS

MULTIFUNCTION CASPAR DISTRACTOR



- Left and right option available
- Locking in 3 planes
- Locking mechanism for Caspar-Pins
- Caspar-Pin threads in different diameters and lengths

COMPATIBLE SYSTEMS

ANTERIOR CERVICAL INTERBODY FUSION (ACIF CAGES)

- SHARK CAGE[®] ACIF
- PINA[®] ECHINUS Ti-coated ACIF
- PINA[®] CERVICUS Invadur[®] ACIF
- PINA[®] CERVICUS Titan ACIF
- PINA[®] VOLUTUS Invadur[®] VBR
- PINA[®] VOLUTUS Titan VBR

Further information at
info@bricon.com
and
+49 (0) 7461 9336-0

BRICON

Bricon[®] GmbH
Eisenbahnstrasse 100
78573 Wurmlingen
Germany

Phone +49 (0) 7461 9336-0

info@bricon.com
www.bricon.com

 0482

bricon.com

Bricon[®] GmbH / All Rights Reserved / Flyer SHARK CP PRO[®] / FL0011_SHARK_CP_PRO_EN_05-2022_Rev1.3