

SHARK CP PRO®

ANTERIOR CERVICAL PLATE






SHARK CP PRO[®]

FEATURES

A great leap on basis of the proven SHARK CP[®] system. Shark CP Pro[®] has a completely new developed locking mechanism ensuring the optimum combination of safety and enhanced usability. The fully optimized instrumentation complement the implant system in setting a new clinical standard for anterior cervical plate.

FEATURES:

-  Innovative locking mechanism
-  Slim waist of 14 mm, maximum width of 16 mm
-  Variable and fixed screws
-  Ultra-low profile with 1,95 mm thickness
-  Screw angulation $\pm 10^\circ$, cranial/caudal up to 25°
-  Large windows for bonegraft application
-  41 Plate sizes, up to 5 Level, length ranging from 12 mm-105 mm







SHARK CP PRO®

IMPLANTS

ARTICLE NUMBER	NEW ARTICLE NO.	DESCRIPTION	LEVEL	LENGTH
700.515.117.1XX	700.515.222.1XX	Anterior Cervical Plate 4	1 Level	12-30mm
700.515.117.1XX	700.515.222.1XX	Anterior Cervical Plate 6	2 Level	30-44 mm
700.515.117.1XX	700.515.222.1XX	Anterior Cervical Plate 8	3 Level	42-66 mm
700.515.117.1XX	700.515.222.1XX	Anterior Cervical Plate 10	4 Level	64-88 mm
700.515.117.1XX	700.515.222.1XX	Anterior Cervical Plate 12	5 Level	75-105 mm
520.200.070		Platecaddy	-	-
520.200.066		Lid for platecaddy	-	-

ARTICLE NUMBER	DESCRIPTION	Ø	LENGTH
700.515.115.41X	Screw, fix, self-drilling	4,0	10-18
700.515.126.11X	Screw, fix, self-tapping	4,0	10-18
700.515.126.12X	Screw, fix, self-tapping, bi-cortical	4,0	20-28
700.515.113.41X	Rescue screw, fix, self-tapping	4,35	11-19
700.515.113.42X	Rescue screw, fix, self-tapping, bi-cortical	4,35	21-29
520.200.068	Screwcaddy, fixed screws	-	-
520.200.069	Lid for screwcaddy	-	-

ARTICLE NUMBER	DESCRIPTION	Ø	LENGTH
700.515.120.11X	Screw, variable, self-drilling	4,0	10-18
700.515.104.11X	Screw, variable, self-tapping	4,0	10-18
700.515.104.12X	Screw, variable, self-tapping, bi-cortical	4,0	20-28
700.515.116.41X	Rescue screw, var., self-tapping	4,35	11-19
700.515.116.42X	Rescue screw, var., self-tapping, bi-cortical	4,35	21-29
520.200.067	Screwcaddy, variable screws	-	-
520.200.069	Lid for screwcaddy	-	-






SHARK CP PRO®

MORE INFORMATION

INSTRUMENTS

MULTIFUNKTIONALER CASPAR DISTRACTOR

- Left and right option available
- Locking in 3 planes
- Locking mechanism for Caspar-Pins
- Caspar-Pin threads in different diameters and lengths



COMPATIBLE SYSTEMS

ACIF CAGES AND VBR

- SHARK CAGE® ACIF
- PINA® ECHINUS Ti-coated ACIF
- PINA® CERVICUS Invadur® ACIF
- PINA® CERVICUS Titan ACIF
- PINA® VOLUTUS Invadur® VBR
- PINA® VOLUTUS Titan VBR

Further Information at
info@bricon.com
and
+49 (0) 7461 9336-0

BRICON

Bricon GmbH
Eisenbahnstrasse 100
78573 Wurmlingen
Germany

☎ +49 7461 9336 - 0
✉ info@bricon.com
🌐 bricon.com



[https://bricon.com/
produkte/wirbelsaeule/](https://bricon.com/produkte/wirbelsaeule/)

